



Single-Use PETG Cap Assemblies

Ace Sanitary provides completely customizable assembly solutions in our Class 7 Cleanroom facility and uses USP Class VI approved materials. Bottles and cap assemblies are also available in multiple brand materials, sizes, connections with no minimum quantities. Clean packed and ready to process with gamma irradiation services for bottles and cap assemblies. Please contact your Ace Sanitary sales representative for more details, additional sizes, or to configure your bottle assembly.

FEATURES

- Cap assemblies come individually double bagged and process ready with labels on the inner bag
- Multiple bottle sizes allow for process-specific flexibility
- Cap assemblies are available in 2 and 3 port configurations with 38-430 neck sizes which can be used on most commercially available bottles
- Cap assemblies in customized size, length, connections, components, with notched dip tubes
- Gamma Irradiation 25kGy to 45kGy for microbial control

APPLICATIONS

Production, Cleanroom, Sampling and Storage, Laboratory, Environmental, Life Sciences

CERTIFICATIONS

USP Class VI, ADFC Animal Derived Component Free, RoHS, USP 87, DEHP Phthalate Free, BPA Free, Latex Free, Non-Pyrogenic, Non-Toxic, DNase/RNase-Free



Multiple end style options are available.

Coupling Expertise with Speed & Simplicity™



Single-Use Cap Assemblies

Part Number	Capacity	Tubing Material	Tubing Length	Cap Size	Connections	Tubing Size
SU/AWT125MLB2PC-GI	125mL, 4oz.	TPE Weldable	Port 1: 18" Port 2: 3"	38-430	Port 1: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/AWT250MLB2PC-GI	205mL, 8oz.	TPE Weldable	Port 1: 18" Port 2: 3"	38-430	Port 1: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/AWT500MLB2PC-GI	500mL, 16oz.	TPE Weldable	Port 1: 18" Port 2: 3"	38-430	Port 1: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/AWT1000MLB2PC-GI	1000mL, 32oz.	TPE Weldable	Port 1: 18" Port 2: 3"	38-430	Port 1: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/AWT125MLB3PC-GI	125mL, 4oz.	TPE Weldable	Port 1: 18" Port 2: 3" Port 3: 18"	38-430	Port 1&3: PP Plug, Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/AWT250MLB3PC-GI	205mL, 8oz.	TPE Weldable	Port 1: 18" Port 2: 3" Port 3: 18"	38-430	Port 1&3: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/AWT500MLB3PC-GI	500mL, 16oz.	TPE Weldable	Port 1: 18" Port 2: 3" Port 3: 18"	38-430	Port 1&3: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/AWT1000MLB3PC-GI	1000mL, 32oz.	TPE Weldable	Port 1: 18" Port 2: 3" Port 3: 18"	38-430	Port 1&3: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/SCT125MLB2PC-GI	125mL, 4oz.	Plt. Cured Silicone	Port 1: 18" Port 2: 3"	38-430	Port 1: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/SCT250MLB2PC-GI	205mL, 8oz.	Plt. Cured Silicone	Port 1: 18" Port 2: 3"	38-430	Port 1: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/SCT500MLB2PC-GI	500mL, 16oz.	Plt. Cured Silicone	Port 1: 18" Port 2: 3"	38-430	Port 1: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/SCT1000MLB2PC-GI	1000mL, 32oz.	Plt. Cured Silicone	Port 1: 18" Port 2: 3"	38-430	Port 1: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/SCT125MLB3PC-GI	125mL, 4oz.	Plt. Cured Silicone	Port 1: 18" Port 2: 3" Port 3: 18"	38-430	Port 1&3: PP Plug, Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/SCT250MLB3PC-GI	205mL, 8oz.	Plt. Cured Silicone	Port 1: 18" Port 2: 3" Port 3: 18"	38-430	Port 1&3: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/SCT500MLB3PC-GI	500mL, 16oz.	Plt. Cured Silicone	Port 1: 18" Port 2: 3" Port 3: 18"	38-430	Port 1&3: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/SCT1000MLB3PC-GI	1000mL, 32oz.	Plt. Cured Silicone	Port 1: 18" Port 2: 3" Port 3: 18"	38-430	Port 1&3: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/APT125MLB2PC-GI	125mL, 4oz.	TPV Biological Pump	Port 1: 18" Port 2: 3"	38-430	Port 1: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/APT250MLB2PC-GI	205mL, 8oz.	TPV Biological Pump	Port 1: 18" Port 2: 3"	38-430	Port 1: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/APT500MLB2PC-GI	500mL, 16oz.	TPV Biological Pump	Port 1: 18" Port 2: 3"	38-430	Port 1: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/APT1000MLB2PC-GI	1000mL, 32oz.	TPV Biological Pump	Port 1: 18" Port 2: 3"	38-430	Port 1: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/APT125MLB3PC-GI	125mL, 4oz.	TPV Biological Pump	Port 1: 18" Port 2: 3" Port 3: 18"	38-430	Port 1&3: PP Plug, Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/APT250MLB3PC-GI	205mL, 8oz.	TPV Biological Pump	Port 1: 18" Port 2: 3" Port 3: 18"	38-430	Port 1&3: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/APT500MLB3PC-GI	500mL, 16oz.	TPV Biological Pump	Port 1: 18" Port 2: 3" Port 3: 18"	38-430	Port 1&3: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"
SU/APT1000MLB3PC-GI	1000mL, 32oz.	TPV Biological Pump	Port 1: 18" Port 2: 3" Port 3: 18"	38-430	Port 1&3: PP Plug Port 2: Acro 37 .2um Filter	3/16"X3/8"

Coupling Expertise with Speed & Simplicity™

